

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Coelenterazine hcp – solid

Product Number: **C3355**
CAS Number: 123437-32-1

Formula: C₂₅H₂₅N₃O₂
Formula Weight: 399.48 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form)	Film
Solubility (Color)	Yellow to Orange
Solubility (Turbidity)	Clear
50 ug + 0.05 mL, MeOH	
Proton NMR spectrum	Conforms to Structure
Purity (HPLC)	≥ 95 %

Specification: PRD.0.ZQ5.10000045837

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.