

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti- β -Actin-Peroxidase antibody produced in mouse - clone AC-15, purified from hybridoma cell culture

Product Number:

A3854**TEST****Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Brown
IgG Concentration (OD 280nm)	2.0 - 4.0 mg/ml
Molar Ratio P/E	0.6 - 1.4
WB-Cell Line/Tissue Extract	Chicken Fibroblasts or HS-68 Human Foreskin Fibroblasts
WB Titer	1:25,000-1:50,000
Recommended Retest Period	-----
1.5 Years	

Specification: PRE.0.ZQ5.10000015199

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.