

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Silver bromide - 99%

Product Number: 226815  
CAS Number: 7785-23-1  
MDL: MFCD00003398  
Formula: AgBr  
Formula Weight: 187.77 g/mol

AgBr

TEST	Specification
------	---------------

Appearance (Color)	Yellow to Yellow-Green
Appearance (Form)	Powder or Crystal or Chunk(s)
Assay (% Ag by NACN) % Ag by Titration With Silver Nitrate After NaCN Treatment	56.5 - 58.3 %
ICP Major Analysis Confirms Silver Component	Confirmed

Specification: PRD.0.ZQ5.10000022184

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.