

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Acetyl- γ -O-hexadecyl-L- α -phosphatidylcholine hydrate - $\geq 98\%$

Product Number:

P4904

CAS Number:

74389-68-7

MDL:

MFCD00150679

Formula:


C₂₆H₅₄NO₇P · xH₂O

Formula Weight:

523.68 g/mol

Storage Temperature:

-20 °C



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear

20 mg/ml, CHCl₃

Purity (TLC)

 $\geq 98\%$

Water (by Karl Fischer)

 $\leq 10\%$

Specification: PRD.1.ZQ5.10000039179

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.