

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Protein kinase C β II isozyme human - >85% (SDS-PAGE), recombinant, expressed in baculovirus infected insect cells, buffered aqueous glycerol solutionProduct Number: **P3287**

MDL: MFCD03458219

Storage Temperature: -70 °C

TEST	Specification
mg protein/ml (Bradford)	0.05 - 1.10
units/mg protein	\geq 500
One unit will transfer 1 nanomole of Phosphate to Histone H3 per min at pH 7.4 at 30 deg C.	
Purity by SDS Electrophoresis	Conforms
> / = 85%	

Specification: PRD.0.ZQ5.10000079770

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.