

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

(11-Mercaptoundecyl)-N,N,N-trimethylammonium bromide

Product Number: 733385

CAS Number: 197587-43-2

Formula: C14H32BrNS Formula Weight: 326.38 g/mol

Storage Temperature: -20 ℃

TEST Specification

Appearance (Color) Appearance (Form) Infrared spectrum Titration by AgNO3

Specification: PRD.1.ZQ5.10000035261

White to Brown Powder or Crystals Conforms to Structure

CH₅ H₅C-N²-CH₂|CH₂|₃CH₅SH Br² CH₅

> 90.0 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.