

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Testosterone - ≥98%

Product Number: T1500 CAS Number: 58-22-0

MDL: MFCD00003654
Formula: C19H28O2
Formula Weight: 288.42 g/mol



TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)
100 mg/mL, CHCl3
1H NMR Spectrum
13C NMR Spectrum
Purity (HPLC)
Recommended Retest Period

3 years

Specification: PRD.1.ZQ5.10000045203

Specification

White to Off White Powder Colorless to Very Faint Yellow

Conforms to Structure Conforms to Structure > 98 %

_

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.