

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Eucalyptus oil - 80-85%, FGProduct Number: **W246611**
CAS Number: 84625-32-1

| TEST | Specification |
|--|--------------------------|
| Appearance (Color) | Colorless to Pale Yellow |
| Appearance (Form) | Liquid |
| Refractive index at 20 ° C | 1.458 - 1.470 |
| Infrared spectrum | Conforms to Structure |
| Specific Gravity | 0.905 - 0.925 |
| Solidification Point | 80.0 - 85.0 |
| (SDPT = 80% to 85% As Cineole C10H18O FCC) | |
| Assay | Pass |
| Passes Phellandrene Test | |
| Arsenic (As) | ≤ 3 ppm |
| Cadmium (Cd) | ≤ 1 ppm |
| Mercury (Hg) | ≤ 1 ppm |
| Lead (Pb) | ≤ 10 ppm |
| Solubility (Turbidity) | Clear |
| 1ML/5ML, 70% ETOH Clear | |
| Clear upon dilution to 10ml | |
| Solubility (Color) | Pass |
| Expiration Date Period | ----- |
| 5 Years | |

Specification: PRD.1.ZQ5.10000062040

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.