

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Agarose - Low EEO

 Product Number:
 A0576

 CAS Number:
 9012-36-6

 MDL:
 MFCD00081294

TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)

1.5 g plus 100 ml H2O

Moisture Content

Residue on ignition (Ash)

Sulfate (SO4)
Gel Strength
1.5% gel
Gel Temperature

1.5% gel Melting Temperature

Turbidity

1.5% gel (Units are expressed in NTU)

Electroendosmosis (-MR)

Recommended Retest Period

3 years

Specification: PRD.2.ZQ5.10000015433

Specification

White to Off White

Pow der

Colorless to Faint Yellow Clear to Very Slightly Hazy

≤ 7 %

< 0.25 %

≤ 0.12 %

 \geq 3200 gm/cm2

34.5 - 37.5 ℃

86.5 - 89.5 ℃

≤ 4

< 0.12

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.