

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Acetonitrile - ReagentPlus®, 99%

Product Number: 110086 CAS Number: 75-05-8

MDL: MFCD00001878

Formula: C2H3N Formula Weight: 41.05 g/mol CH₃CN

TEST Specification

Appearance (Color) Colorless

Appearance (Form) Liquid

Infrared spectrum Conforms to Structure

Purity (GC) \geq 98.5 %

Specification: PRD.1.ZQ5.10000024662

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.