

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium taurohyodeoxycholate hydrate - ≥98%

Product Number: T0682
CAS Number: 38411-85-7
MDL: MFCD09263185
Formula: C₂₆H₄₄NNaO₆S·xH₂O
Formula Weight: 521.69 g/mol



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	
Water (by Karl Fischer)	≤ 7 %
Purity (TLC)	≥ 98 %
Carbon (anhydrous)	54.5 - 60.0 %
Nitrogen (anhydrous)	2.0 - 3.4 %

Specification: PRD.0.ZQ5.10000027161

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.