


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Glycerol – ReagentPlus[®], ≥99.0% (GC)

Product Number: G7757
CAS Number: 56-81-5
MDL: MFCD00004722
Formula: C₃H₈O₃
Formula Weight: 92.09 g/mol



TEST	Specification
Appearance (Clarity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Viscous Liquid
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 99.0 %
Water (by Karl Fischer)	≤ 1.0 %
Heavy Metals (as Lead)	≤ 5 ppm
Sulfated Ash	≤ 0.1 %
Recommended Retest Period	-----
3 Years	

Specification: PRD.0.ZQ5.10000020931

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.