

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-ATP Synthase β antibody produced in mouse - 1 mg/mL, clone 4.3E8.D10, purified immunoglobulin, buffered aqueous solutionProduct Number: **A9728**

MDL: MFCD03455681

TEST	Specification
------	---------------

Formulation	Conforms
Solution in phosphate buffered saline containing 1 mg/mL BSA and 0.05% sodium azide	
Concentration	Conforms
1 mg/mL	
Working Dilution	Conforms
By Immunofluorescence; 1:100 - 1:1000	

Specification: PRD.0.ZQ5.10000035585

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.