

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Collagen Type IV from human cell culture – Bornstein and Traub Type IV, 0.3 mg/mL, sterile-filtered,
BioReagent, suitable for cell culture

Product Number: C6745

CAS Number: 9007-34-5

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Solution
Sterility	Pass
Protein Content	0.2 - 0.4 mg/ml
In 25% Acetic acid	
Purity by SDS Page	Conforms
Collagen type IV bands detected between 170-180 KDA	
Country of Origin	United States of America
Country of Manufacture	United States of America
Formulation	Conforms
in 25% HOAc	
Tested For Infectious Agents	Tested

Specification: PRD.1.ZQ5.10000013288

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.