

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lead(II) oxide - 99.999% trace metals basis

Product Number: 203610
CAS Number: 1317-36-8
MDL: MFCD00011164
Formula: OPb
Formula Weight: 223.2 g/mol

PbO

TEST	Specification
-------------	----------------------

Appearance (Color)	Yellow to Orange
Appearance (Form)	Powder or Powder with Chunk(s)
Titration with EDTA % Pb	91.4 - 94.2 %
ICP Major Analysis Confirms Lead Component	Confirmed
Purity 99.999% Purity Based On Trace Metals Analysis	Conforms
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.1.ZQ5.10000066981

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.