

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Collagen-Fluorescein bovine - collagenase substrate

Product Number: C5858
MDL: MFCD02264428
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Faint Yellow to Dark Yellow
Appearance (Form)	Fibers or Solid
Suitability	Pass
Suitable as fluorescent substrate for collagenase	
Specification: PRD.1.ZQ5.10000034869	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.