

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Interleukin-7 Receptor  $\alpha$ /CD127 antibody produced in goat - affinity isolated antibody, lyophilized powderProduct Number: **I8409**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Lyophilized Powder
Western Blot	Conforms
0.1 ug/ml recommended concentration	
Flow Cytometry	Conforms
2.5 $\mu$ g/10 <sup>6</sup> recommended concentration	

Specification: PRD.0.ZQ5.10000083288

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.