

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Ammonium acetate solution – Reagents designed and manufactured under cGMP controls suitable for use in IVD application.

Product Number:

SRE0035



Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Bioburden	≤ 100 CFU/ml
DNase, Exonuclease Detection	None Detected
RNAse Detection	None Detected
NICKase, Endonuclease Detection	None Detected
Protease Detection	None Detected
Alkaline Phosphatase Impurity	None Detected
Expiration Date Period	Confirmed
2 years	

Specification: PRD.0.ZQ5.10000101271

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.