

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium-potassium alloy - Potassium 78 wt. %, Sodium 22 wt. %

Product Number:


244872

CAS Number:

11135-81-2

MDL:

MFCD00134088



TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Silvery Metallic

Appearance (Form)

Liquid

Potassium (K)

76.0 - 80.0 %

Trace Metal Analysis

≤ 550 ppm

Balance Sodium

Purity

Meets Requirements

99.95% Based on Trace Metals Analysis

Specification: PRD.1.ZQ5.10000008234

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.