


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Silicone oil – viscosity 500 cSt (25 °C)

Product Number: **378380**
CAS Number: 63148-62-9
MDL: MFCD00132673
Formula: $(-\text{Si}(\text{CH}_3)_2\text{O}-)_n$



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Viscous Liquid
Infrared spectrum	Conforms to Structure
Viscosity At 25°C	470.0 - 530.0 cps
Specification: PRD.2.ZQ5.10000005725	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.