


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Vincamine - 98%

Product Number: **194719**
CAS Number: 1617-90-9
MDL: MFCD00078054
Formula: C₂₁H₂₆N₂O₃
Formula Weight: 354.44 g/mol

**TEST****Specification**

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder or Crystals or Chunk(s)
Infrared Spectrum	Conforms to Structure
Titration with HClO ₄	97.5 - 102.5 %
Optical Rotation c = 1%, Pyridine	39.8 - 45.8 %

Specification: PRD.1.ZQ5.10000043420

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.