

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Barium acetate - 99.999% trace metals basis

Product Number:

255912

CAS Number:

543-80-6

MDL:


MFCD00012447

Formula:

C₄H₆BaO₄

Formula Weight:

255.42 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder or Crystal or Chunk(s)

ICP Major Analysis

Confirmed

Confirms Barium Component

Purity

Meets Requirements

99.999% Based On Trace Metals Analysis

Trace Metal Analysis

≤ 15.0 ppm

Specification: PRD.1.ZQ5.10000034035

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.