

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,2-Didecanoyl-sn-glycero-3-phosphocholine ≥99% (TLC)

Product Number:

P7081

CAS Number:

3436-44-0

MDL:

MFCD00063538

Formula:


C₂₈H₅₆NO₈P

Formula Weight:

565.72 g/mol

Storage Temperature:

-20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Pow der
Solubility (Color)	Colorless
Solubility (Turbidity) 25 mg in 0.5 mL, CHCl ₃	Clear
Proton NMR spectrum	Conforms to Structure
Purity (TLC)	≥99 %
Purity (GC) As Decanoic Acid	≥99 %
Recommended Retest Period 5 Years	-----

Remarks:

Specification Date : 02/10/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.