

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Castor oil – meets testing specifications per USP

Product Number:

C9606

CAS Number:

8001-79-4

MDL:

MFCD00130746

TEST	Specification
Specific Gravity	0.957 - 0.961
Distinction	Pass
From Most Other Fixed Oils	
Heavy Metal	≤ 0.001 %
Free Fatty Acids	Pass
Hydroxyl Value	160.0 - 168.0
Iodine Value	83 - 88
Saponification Value	176.0 - 182.0
Meets USP Test Specifications	Current USP Specifications

Specification: PRD.1.ZQ5.10000034551

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.