

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Palladium(II) chloride - ReagentPlus[®], 99%

Product Number: 205885
CAS Number: 7647-10-1
MDL: MFCD00003558
Formula: Cl₂Pd
Formula Weight: 177.33 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Brown to Very Dark Brown and Brown-Red and Red-Brown and Very Dark Red-Brown
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Palladium (Pd)	59.1 - 60.9 %
Solubility (Turbidity)	Clear
1% In 5 N HCL	
Solubility (Color)	Brown to Dark Brown and Brown-Red to Dark Brown-Red

Specification: PRD.1.ZQ5.10000035461

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.