SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: (2-Hydroxypropyl)-β-cyclodextrin - average Mw ~1,380

Product Number:	332593
CAS Number:	128446-35-5
MDL:	MFCD00069372



TEST	Specification	
Appearance (Color)	White	
Appearance (Form)	Powder	
Infrared spectrum	Conforms to Structure	
Ratio	0.5 - 0.7	
Molar Substitution Ratio Average		
Optical Rotation (Degrees)	124 - 154	
Degrees		
C= 1%, H2O		

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.