

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-McCP2 antibody produced in rabbit - ~0.6 mg/mL, affinity isolated antibody, buffered aqueous solution

Product Number:

**M9317****TEST****Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	0.6 - 0.8 mg/ml
WB-Cell Line/Tissue Extract	MCF7 Hu Breast Carcinoma
WB Working Concentration	0.5 - 1.0 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000017308

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.