

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Yeast Carbon Base

Product Number: Y3627

Storage Temperature: 2 - 8 ℃

TEST Specification

Appearance (Color) White to Off-White

Appearance (Form) Powder

Homogenous, free-flowing

pH 5.30 - 5.70

1.17% w/v aqueous solution at 25°C

Solubility (Turbidity) Clear
Solubility (Color) Colorless

Solubility (Form) Solution

1.17% w/v aqueous solution at 25°C

Cultural Characteristics Pass

Pass as described in data sheet.

Cultural characteristics observed after an incubation at 25-30°C for 6-7 days or longer if required

Saccharomyces cerevisiae Conforms

Growth(plain): none - poor

Growth (with Ammonium Sulfate): good

Saccharomyces uvarum Conforms

Growth (plain): none - poor

Growth (with Ammonium sulfate): good

Specification Date: 03/30/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.