


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Palmitoyl-2-oleoyl-3-linoleoyl-rac-glycerol - ≥98%

Product Number: P8577
CAS Number: 2680-59-3
MDL: MFCD00214386
Formula: C55H100O6
Formula Weight: 857.38 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
% Purity (TLC)	≥ 98

Specification: PRD.1.ZQ5.10000061558

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.