SIGMA-ALDRICH®

sigma-aldrich.com


3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Gold nanoparticles - 50 nm diameter, methyl terminated, PEG 2000 coated, OD 50, dispersion in H2 O ()

Product Number:

765678



Storage Temperature: 2 - 8 °C

Dark Red

Specification: PRD.0.ZQ5.10000027704

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.