

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Citraconic anhydride 98%

Product Number: 125318
CAS Number: 616-02-4
MDL: MFCD00005522
Formula: C₅H₄O₃
Formula Weight: 112.08 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥97.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
Recommended Retest Period	-----
2 Years	

Remarks:

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.