

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-GABAA Receptor (δ subunit), N-Terminal antibody produced in rabbit – affinity isolated antibody, buffered aqueous glycerol solution

Product Number:

G9794

MDL:

MFCD07370325

TEST	Specification
------	---------------

Formulation	Conforms
Solution in 10 mM HEPES, pH 7.5, containing 150 mM NaCl, 100 ug/mL BSA and 50% Glycerol	
Volume	Conforms
100 uL	
Working Dilution by Direct Imm	Conforms
1:1,000 (recommended dilution)	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.