

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
6-Mercapto-1-hexanol - 97%

Product Number: 451088
CAS Number: 1633-78-9
MDL: MFCD00068552
Formula: C₆H₁₄OS
Formula Weight: 134.24 g/mol

SHCH₂(CH₂)₄CH₂OH

TEST	Specification
------	---------------

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Titration with KIO₃

96.5 - 103.5 %

Purity (GC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000014062

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.