

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silicone oil - for melting point and boiling point apparatuses

Product Number:

146153


CAS Number:

63148-62-9

MDL:

MFCD00132673

Formula:

 $(-\text{Si}(\text{CH}_3)_2\text{O}-)_n$ **TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Viscous Liquid

Refractive index at 20 ° C

1.402 - 1.406

Infrared spectrum

Conforms to Structure

Viscosity

45.0 - 55.0

At 25°C; Centistokes

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.