

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Chloropyrimidine - 95%

Product Number: 193291
CAS Number: 1722-12-9
MDL: MFCD00006060
Formula: C₄H₃ClN₂
Formula Weight: 114.53 g/mol

TEST	Specification
Appearance (Color) White to Yellow or Orange	Conforms to Requirements
Appearance (Form)	Crystals
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 94.5 %
Note May Contain 2-Hydroxypyrimidine	Conforms

Specification: PRD.2.ZQ5.10000004390

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.