


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tantalum(V) ethoxide - packaged for use in deposition systems

Product Number: **760404**
CAS Number: 6074-84-6



Formula: C₁₀H₂₅O₅Ta
Formula Weight: 406.25 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow and Colorless to Orange
Appearance (Form)	Liquid
Proton NMR Spectrum	Conforms to Structure
Exists as a Bridged Dimer	

Specification: PRD.2.ZQ5.10000021099

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.