

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


L- α -Phosphatidyl-DL-glycerol sodium salt from egg yolk lecithin - $\geq 99\%$ (TLC), lyophilized powder

Product Number:

P8318

MDL:

MFCD00213550



Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Very Faint Yellow

Solubility (Turbidity)

Clear

25 mg/mL, CHCl₃

Phosphorus (P)

3.5 - 4.5 %

Sodium (Na)

1.6 - 3.2 %

Purity (TLC)

 ≥ 99 %

Specification: PRD.6.ZQ5.10000006867

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.