

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

UbcH9 Conjugating Enzyme human - ≥95% (SDS-PAGE), recombinant, expressed in E. coli, buffered aqueous

glycerol solution

Product Number: U9257

MDL: MFCD03458912

Storage Temperature: -70 ℃

TEST	Specification	
Formulation Solution in 50 mM HEPES, pH 7.5, 100 mM NaCl, 10 % glycerol, 1 mM TCEP		
Protein Content	≥ 0.5 mg/ml	
Purity by SDS Electrophoresis	≥ 95 %	
Note		

Specification: PRD.1.ZQ5.10000048377

BIOCHEM, INC

MANUFACTURED FOR SIGMA BY BOSTON

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.