

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chromium(II) chloride - 95%

Product Number: 244805
CAS Number: 10049-05-5
MDL: MFCD00010947
Formula: Cl₂Cr
Formula Weight: 122.90 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Light Green to Green-Grey	
Appearance (Form)	Powder
Titration with AgNO ₃	54.5 - 60.9 %
% Cl	
ICP Major Analysis	Confirmed
Confirms Chromium Component	

Specification: PRD.1.ZQ5.10000025845

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.