


Product Specification

Product Name:
Clarithromycin - $\geq 95\%$ (HPLC)

Product Number: **C9742**
CAS Number: 81103-11-9
Formula: C₃₈H₆₉NO₁₃



TEST	Specification
Appearance (Form)	Powder
Appearance (Color)	Conforms
White to off white	
Purity (HPLC)	$\geq 95\%$
Identification IR	Conforms
Recommended Retest Period	-----
4 Years	

Specification: PRE.5.ZQ5.10000021227

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

