

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Vinyl decanoate – contains monomethyl ether hydroquinone as inhibitor, >99%

Product Number:

411795

CAS Number:

4704-31-8

MDL:


MFCD00027349

Formula:

C₁₂H₂₂O₂

Formula Weight:

198.30 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 ° C	1.434 - 1.438
Purity (GC)	≥ 98.5 %
Inhibitor	Confirmed
Monomethyl ether hydroquinone	

Specification: PRD.1.ZQ5.10000010277

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.