

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

6-Benzyloxy-5-methoxy-2-carboxyindole - ≥99%, powder

Product Number:

B1251

CAS Number:

2495-92-3

MDL:

MFCD00055977

Formula:


C₁₇H₁₅NO₄

Formula Weight:

297.31 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

Off-White to Yellow

Appearance (Form)

Powder

Proton NMR spectrum

Conforms to Structure

Carbon (anhydrous)

67.9 - 69.3 %

Nitrogen (anhydrous)

4.6 - 4.8 %

Water (by Karl Fischer)

≤ 3 %

Purity (TLC)

≥ 99 %

Solubility (Color)

Yellow to Orange

Solubility (Turbidity)

Clear

50 mg/ml, pyridine

Specification: PRD.0.ZQ5.10000052017

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.