

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
PIPES - ≥99% (titration)

Product Number: P6757
CAS Number: 5625-37-6
MDL: MFCD00006159
Formula: C₈H₁₈N₂O₆S₂
Formula Weight: 302.37 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Crystalline Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
20 g plus 80 mL, 1 N NaOH	
Infrared spectrum	Conforms to Structure
UV Absorbance 290nm	≤ 0.1
20% W/W in 1 N NaOH	
ICP-MS	≤ 5 ppm
(Total: Ag, As, Bi, Cd, Cu, Hg, Mo, Pb, Sb, Sn)	
Titration with H ₂ SO ₄	≥ 99 %
Recommended Retest Period	-----
3 years	

Specification: PRD.3.ZQ5.10000038534

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.