

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Nicotinamide adenine dinucleotide, reduced disodium salt - ~98%, pkg of 25-30 mg (per vial)

Product Number:

N9535

CAS Number:

606-68-8

MDL:


MFCD00036200

Formula:

C₂₁H₂₇N₇Na₂O₁₄P₂

Formula Weight:

709.4 g/mol

**TEST****Specification**

Appearance (Color)

White to Yellow

Appearance (Form)

Powder

Solubility (Color)

Colorless to Yellow

Solubility (Turbidity)

Clear

Solubility (Form)

Solution

mg/vial

25 - 30

Suitability

Suitable

Suitable for use as a co-factor in enzyme assays. Gives comparable response as previous lots in LDH enzyme assay system.

Specification: PRD.2.ZQ5.1000009348

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.