


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Pyrilamine maleate salt

Product Number: **P5514**
CAS Number: 59-33-6
MDL: MFCD00069333
Formula: $C_{17}H_{23}N_3O \cdot C_4H_4O_4$
Formula Weight: 401.46 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Titration by HClO ₄	≥ 98 %
Purity (TLC)	≥ 98 %
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
100 mg/mL, H ₂ O	
Recommended Retest Period	-----
4 years	

Specification: PRD.0.ZQ5.10000074555

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.