

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Peru balsam oilProduct Number: **W211710**
CAS Number: 8007-00-9
MDL: MFCD00147893

TEST	Specification
Appearance (Color)	Pale Yellow to Brown
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.565 - 1.575
Infrared spectrum	Conforms to Structure
Purity	Pass
%, Similar to Standard	
Specific Gravity	1.110 - 1.120
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.0.ZQ5.10000037985

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.