

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Lymphotoxin β Receptor antibody produced in goat - affinity isolated antibody, lyophilized powderProduct Number: **L5412**

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Lyophilized from 0.2 um-filtered solution in phosphate buffered saline containing carbohydrates	
Endotoxin Assay	Conforms
< 0.1 EU/ug antibody	
Working Dilution	Conforms
By Indirect Immunoblotting; 0.1 - 0.2 ug/mL using human lymphotoxin beta R at 5 ng/lane under non-reducing and 25 ng/lane under reducing conditions.	
Flow Cytometry	Conforms
2.5 ug antibody/10 ⁶ cells	

Specification: PRD.0.ZQ5.10000070478

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.