


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-(Benzyloxy)aniline hydrochloride - 98%

Product Number: **116637**
CAS Number: 51388-20-6
MDL: MFCD00012995
Formula: C₁₃H₁₃NO · HCl
Formula Weight: 235.71 g/mol



TEST

Specification

| | |
|--------------------------------|---|
| Appearance (Color) | Off White to Beige and Faint Brown to Light Brown |
| Appearance (Form) | Powder |
| Proton NMR spectra | Conforms to Structure |
| Titration by AgNO ₃ | 97.5 - 102.5 % |

Specification: PRD.0.ZQ5.10000041033

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.