

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Albumin from chicken egg white - lyophilized powder, $\geq 98\%$ (agarose gel electrophoresis)

Product Number: A2512
CAS Number: 9006-59-1
MDL: MFCD00130424

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/ml, H ₂ O	
Water (by Karl Fischer)	$\leq 5\%$
Nitrogen	14.0 - 16.0 %
mol/mol Albumin	4 - 20
Mannose content	
Agarose Electrophoresis	$\geq 98\%$
Recommended Retest Period	-----
2 years	

Specification: PRD.2.ZQ5.10000028815

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.