

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Nuclease S1 from *Aspergillus oryzae* – for single-strand DNA/RNA digestion

Product Number:

N5661

CAS Number:

37288-25-8

MDL:

MFCD00131010

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/ml

≥ 100000

One unit will cause 1.0 microgram of single-stranded nucleic acid to become perchloric acid soluble per minute at pH 4.6 at 37 deg C.

Specification: PRD.0.ZQ5.10000045272

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.